Philips Speech Mike Driver Extension Setup

Citrix XenApp does not, by default, pass through the record button on the **Philips Speech Mike III**. However, the button press will work with the manufactures driver extension installed.

Installing

Note: For purposes of our documentation, we are using Windows 10. Windows 7 or higher is supported with this version of the Philips Speech Mike Driver extensions.

Download the Philips Speech Mike Driver installation package from the following location:

http://support.pinnaclehealth.org/PhilipsSpeechDriversSetup12-3.exe

Right Mouse Click the Installer. Select Run as Administrator.



The InstallShield Wizard Starts.



Ensure ONLY the Client Extension for Citrix is selected.



Click Next.

Ignore DPM Mounting Parameters. Click Next.

Philips Speech Drivers - InstallShield Wizard	×
DPM mounting parameters Please specify DPM mounting parameters	
In DPM Local path, please enter a valid folder on a local hard disk of this computer. This path will show the contents of the DPM memory card. This path will be mapped to the server. DPM Local Path	
CADPM In DPM Server path please, ever the path where the DPM Local path is mapped by Citric. If Drive remapping is enabled on the Citrix server then the two paths will be the same, otherwise the C: drive on the client will become V: on the server and D: will become U: and so on. Please, make sure to enter the appropriate value	
DPM Server Path C:\DPM	
InstellShield	

Click Install then Click Finish.

Note: A reboot of the Windows Endpoint is not required if the Citrix Receiver was not running at the time of installation. A reboot is recommended prior to testing if possible.

Note: You may see an HDX Warning that a Virtual Desktop is attempting to use your microphone or webcam. **Click** Permit use of these devices.